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High-Quality Protein Crystallization Project for Protein Structure and Function Analyses in Support of Practical Applications

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Protein Crystallization Experiments in Space


To prepare for the full-scale use of Kibo, the Japanese Experiment Module (JEM) that will be attached to the International Space Station (ISS), the Japan Aerospace Exploration Agency (JAXA) will perform a series of six experiments in space from 2003 through 2005. These experiments grow protein crystals aboard the Service Module of the ISS.

JAXA is conducting this project in cooperation with research organizations and universities that analyze the structures and functions of proteins as part of a national project. It is also working with private-sector companies that seek to create genome-based pharmaceuticals. It is difficult to grow good single crystals of proteins on the ground that are suitable for more precise analyses. In such cases, crystal growth experiments will be conducted in a space environment. Afterward, the protein structures and functions will be analyzed on the ground, and the results will be utilized in practical applications.

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Relationship between space and the function and structure analyses of protein crystals

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Significance of Protein Structure and Function Analyses

In April 2003, it was announced that the human genome had been completely decoded. Research has now moved to the next step, post-genome research, which involves identifying the structures and functions of proteins made based on genes. In order to learn more about protein functions, it is extremely important to clarify the three-dimensional protein structure at the atomic level. In most cases, the structure of a protein is closely related to its function, so in-depth analysis of the three-dimensional protein structure is being linked to the study of protein functions.

Human life is supported by nearly 100,000 different kinds of proteins, most of whose structures are still unknown. By analyzing the structures of many different kinds of proteins and understanding how they interact with one another, it may become possible to identify various mechanisms involved in life phenomena. It should also become possible to inhibit or activate functions of certain proteins to develop new pharmaceuticals.

Pharmaceutical companies are becoming increasingly interested in using the results from post-genome research to develop new medicines using Structure-Based Drug Design (SBDD).

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
Through the first, second and third space experiments, JAXA launched proteins of alpha-amylase and lysozyme to validate the developed crystallization technique. As a result, high-quality protein crystals with the best resolutions in the world were successfully obtained, and the structure was analyzed closely and carefully.

In some of the proteins (e.g., sleeping and allergy material synthetase, and proteins related to symptoms of parasite infection) provided by user organizations, we also obtained the highest resolution data from the fine protein crystal grown in space; this is expected to lead to pharmaceutical development in the future.

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
ISS is an "in-space laboratory" composed of many modules. A variety of space experiments are currently performed in Russian Service Module and US Experiment Module.

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
The Russian Service Module "Zvezda" (meaning "star" in Russian) where GCF experiments are conducted was launched in July 2000. Zvezda provides essential functions aboard the ISS, such as life support, flight control, and communications with flight controllers on the ground. It also enables automatic docking to Soyuz spaceships and Progress supply ships.

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The Progress supply ship (an improved Soyuz spaceship) is an automated, unmanned transport ship. It transports GCF as well as fuel, oxygen, water, food, clothing, spare units, and mail. Progress functions as a propulsion system for the ISS by boosting its orbit, after docking to the ISS.

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Soyuz is a manned Russian spacecraft carrying crew members and cargo from/to the ISS. One Soyuz is docked to the ISS at all times to serve as an emergency escape vehicle. A new Soyuz is launched every six months and replaces the old Soyuz. Once the on-orbit experiment is completed, a Soyuz returns the GCF to the ground.

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Protein Crystallization in a Space Environment

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
X-ray crystallographic analysis, in which X-rays are reflected and scattered from a protein crystal, is commonly used to examine the three-dimensional (3-D) structure of large-molecule proteins. X-ray diffraction patterns obtained from the reflection and scattering are protein-unique, so the 3-D structure of the protein crystal can be estimated by analyzing such patterns. More accurate estimation requires that crystals have molecules arrayed in good order. The better the molecules are arrayed, the clearer X-ray diffraction patterns become, and therefore, the more precisely the three-dimensional structure can be estimated.


It is difficult to produce high-quality protein crystals on the ground due to convection and sedimentation caused by gravity, which disturb the molecule array by causing crystals to stick to each other. In the space environment, however, high-quality crystals with molecules arrayed in good order can be grown in the near-zero gravity environment that eliminates convection and sedimentation. JAXA has succeeded in growing high-quality crystals, throughout the on-orbit protein crystallization experiments of past missions.

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
Experimental Facility for Growing High-Quality Protein Crystals

The High-Quality Protein Crystallization Project utilizes the Granada Crystallization Facility (GCF) developed jointly by the European Space Agency (ESA) and Granada University in Spain, and the JAXA Crystallization Facility (JCF) developed by JAXA. GCF and JCF contain 23 Granada Crystallization Boxes (GCBs) at maximum and a temperature data logger.

Proteins are crystallized in space inside an incubator on board the Service Module. When samples are returned to the ground by Soyuz, insulators wrapped around the GCF and JCF, a thermos-like container, isolate the crystallization inside GCBs from fluctuations of ambient temperature.

One GCB has up to 11 capillaries, so 11 different protein crystallization experiments at maximum can be conducted at the same time.

結晶生成セル
(サイズ:36×101×7mm)
Granada Crystallization Box
(Dimensions:33×100×7mm)

JAXA結晶生成装置
(サイズ:φ90×165mm)
JAXA Crystallization Facility
(Dimensions:φ90×165mm)


グラナダ結晶生成装置
(サイズ:134×134×84mm)
Granada Crystallization Facility
(Dimensions:134×134×84mm)

真空断熱材 (GCF用)
Vacuum Insulator (for GCF)

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Crystallization Method used in GCBs (Counter diffusion method)

GCB adopts the counter diffusion method that promotes diffusion of protein solutions and crystallization agent solutions in opposite directions inside capillaries through a gel layer. Protein solution diffuses outward; agent solution diffuses inward. This forms a slope of the two concentrations, and crystallization starts at the point and time where the two concentrations meet the optimum condition for crystallization. The tube, filled with gel and fitted to the bottom of the capillary, moderates the diffusion velocity.

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
Overview of Experiments

To promote this project, JAXA is forming cooperative agreements and signing joint contracts with organizations and universities researching the structures and functions of proteins as part of a national project. JAXA is also seeking cooperation with private-sector companies striving to create genome-based pharmaceuticals.

First, protein crystallization conditions are investigated using the counter-diffusion method based on ground conditions provided by the contracted research groups. Appropriate experiments are then conducted on the ground. The results are subsequently used to determine the conditions for filling the GCBs with protein and crystallization agent solutions for the experiments in space. After that, the GCBs are loaded onto the GCF.

With the support of RSC Energia Co., Ltd., of Russia, this experiment facility is launched into space aboard a Progress supply ship from the Baikonur Cosmodrome in Kazakhstan. When the GCF reaches the ISS, it is installed inside the Service Module where it is used to grow protein crystals for two to four months. The GCF is then loaded onto a Soyuz bound for Kazakhstan.

The protein crystals in the GCF are shipped to Japan, where they are divided among various contracted research groups. They obtain data needed for structural analyses, assess the quality of the crystals, and analyze the three-dimensional structure of the protein.

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
第1回目の軌道上実験で得たタンパク質結晶
Protein crystal obtained during the flight of JAXA's first crystallization experiment

Future Plans

JAXA seeks to make significant contributions to studies in the fields of protein structure and function analyses and pharmaceutical development by providing more opportunities to various user organizations such as private companies and national research institutes and universities.

To enable as many users as possible to participate in-space crystallization experiments, JAXA will continue to promote experiment facility development and to verify the in-space crystallization technique. It will also develop support systems that flexibly respond to users' needs.

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X-ray diffraction image

Compared with crystals made on the ground, diffraction images of high-quality crystals grown in space have more as well as clearer diffraction points. This can provide more detailed data closely linked with three-dimensional structural analyses.

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Home Page URL
<http://www.jaxa.jp/>
<http://iss.sfo.jaxa.jp/>

Japan Aerospace Exploration Agency
Office of Space Flight and Operations
1-1, Sengen 2-chome, Tsukuba-shi Ibaraki 305-8505, Japan
Phone:#81-29-868-3074 Fax:#81-29-868-3950

2004.11.10現在
Updated on Nov. 10, 2004